




Laparoscopy

DCS-103E Pro


All in one system for Laparoscopy, Endoscopy, Proctology




Cholecystectomy



Colorectal Cancer




Pancreatic Cancer




Hepatectomy

LAPAROSCOPES 10 MM



LAPAROSCOPES 5 MM



Laparoscope



CCD Coupler



f = 18.5mm

Optimal solution for Telemedicine system



Full HD(1080i / 60 Frames)

Convenient joystick




Image Conversion





Image management



※ Zoom IN/OUT (X1.0 - X3.0)

Easy Connection with PC/Laptops



Functions	FULL HD Live Image (16:9) Image Comparison USB2.0 PC Slave Mode Image Adjustment
Camera Resolution	1920x1080
Video Output	1920x1080p, 60 frame, HDMI
PC Interface	USB 2.0 / Slave Mode Directly connected to the PC application
Internal Storage	32GB
Enternal Storage interface	USB 2.0 Memory Stick / FAT32
Still image Capture	1920x1080 jpg File
User Key Interface	5 Directional Key
Light Source	18W LED lamp, Color temperature (4800. K), Adjustable light intensity
Cable Length	2.0M
Foot switch	for Freeze & Still Capture
AC Input	90~264VAC Max 45Watt Universal AC Input / Full range
Dimension	Width(209mm) x Depth(101mm) x Height(250mm)



AMT, Plot No- A815, MIDC Industrial Area,
Kopar Khairane, Navi Mumbai, Maharashtra 400705

Office : +91-9322673077 | +91-9892703077

Office : +91-8080321777 | +91-9321548914

Web : www.amtmedicals.com

Email : info@amtmedicals.com