



COSMETIC GYNECOLOGY

Equipment




BENEFIT :

- Minimally Invasive, No Incision
(no cutting, no bleeding, no ablation)
 - Non-anesthesia, Painless
 - Don't need any special care before & after procedure
 - No Consumables
 - Easy and Fast procedure
- High rate of success with only one treatment

SPECIFICATIONS :

Laser Type	: CO ₂ laser
Wave Length	: 10.6 um
Watt	: 1 - 50W
Electric Requirements	: 220 - 240 VAC, 50/60Hz
Dimension	: 356(W) X 428(D) X 1260(H) mm
Weight	: 45 Kg



LASER VAGINAL TIGHTENING (LVT)

DIFFERENT MODES :

- | | |
|-------------------|---|
| Surgical Mode | : Surgery, Nevus, Wart, Tumor removal |
| Zoom Mode | : Similar to surgical Application, with bigger spot size |
| Fractional Mode | : Aesthetic, Acne scar, Pore, Wrinkle, Laser resurfacing, Pudenda whitening
Skin/Pudenda texture improvement |
| Beauty Pulse Mode | : Vaginal rejuvenation, Vaginal elasticity Improvement, Colpoxerosis and SUI improvement |



Surgical



Fractional




Vaginal



Zoom

1470nm / 15 W
Diode Laser with LVR and SUI Probe
For Cosmetic Gynecology



SUI (Stress Urinary Incontinence) Handpiece

THE BENEFITS OF LASER TREATMENT

A Simple Non Surgical Process that result in

- Increased Collagen
- Improved Muscle Tone
- Improved Sensation
- Reduced Vaginal Dryness
- Enhanced Orgasms
- OPD process without anesthesia
- No side effects
- Effective and painless

APPLICATIONS

- Vaginal rejuvenation
- GSM-Genitourinary syndrome of menopause
- Vaginal atrophy
- Vaginal laxity
- Diseases associated to postpartum alteration
- HPV
- Dryness
- Itching
- Vaginismus

SPECIFICATION:

● Laser type	: Diode, Semiconductor	● Power Supply of the Laser	: DC 24V/8.33A from the Separate AC
● Wavelength	: 1470nm	● Power Supply of AC Adapter	: Single phase 100-240 VAC: 50-60 Hz;
● MAX Power	: 15 Watts	● AC Adapter	: DC 24v/8,33A Medical device class
● Aiming Beam	: 635nm/Pulse, Single, Group		EN60601-1 Class
● Operating Mode	: Continuous Pulse, Single, Group	● Laser Dimensions	: 27cm x 24,5cm x 9cm
● Pulse Time	: 0.05ms-1000ms	● Laser weight	: 2.5kg.
● Beam delivery	: Fiber Optic with SMA 905 Connector	● Laser Case Dimensions	: 53cm x 38cm x 23cm
● Optical Fiber Connector	: SMA 905, accepts optical fibers having a core diameter from 200µm up to 600µm,	● Weight of Laser with Cases	: Approx 10 kg
		● Class of Medical Device	: IIB
● Emmition Initiation.	: Footswitch	● Laser Safety Class	: 4
● Controller	: Microprocessor	● Electric Safety Class	: I type B
● Display	: 8 TFT with touch panel	● Environmental Conditions During Work	: from +10°C to +25°C Relative humidity
● Cooling System	: Internal, air and thermoelectric cooling		from 30% up to 60%

RF (Radio Frequency)

- Vaginal Rejuvenation
- Majora/ Minora Development
- Skin Lustering



SPECIFICATION:

Input Power	: 100-120/200-240V \pm 10%, 50/60 Hz
Output Power	: 20 watts (maximum) 50-300 Ω 0-10 volts, in Stimulate Mod
Maximum Output Voltage	: 80V, RMS
Operating Frequency	: 460 kHz (\pm 5kHz)
Set Temperature Range	: 35-47°C; this range varies for each Power Delivery Mode
Dimensions	: 5.75" H X 12.5" W X 13.0" D (146mm X 318mm X 330mm)
Weight	: 12.6 lbs (5.7kg)

HIFU (HIGH INTENSITY FOCUSED ULTRASOUND)



- 1) Sonic Vaginal rejuvenation
- 2) Sonic Vaginal tightening
- 2) Vaginal pressure measurement

PLASMA (Centrifuge + Plasma Generator + Needle)



Centrifuge



Bio Filler

PRP Kit



THIS BOX CONTAINS

- 3No - 15ml Tubes
- 1No - Pipette
- 1No - Ampule
- 1No - Blood Removal Cannula
- 1No - Syringe (3 ml) with Needle





SPECIFICATION:

Laser Type	: Diode Laser
Display System	: 10.4" 24 Color touch screen
Wavelength	: 808nm/755+808+1064nm 3 waves
Laser handle shots	: 50,000,000 unit
Energy Scope	: 1-80J/cm (Handle 300W) 1-139J/cm (Handle 500W/600W) 1-223J/cm (Handle 800W) 1-278J/cm (Handle 1200W)
Pulse Duration	: 1-400ms adjustable
Frequency	: 1-15 HZ
Cooling system	: Tec Condenser + Air + Water
Machine Power	: 3000W
Handle Power	: 3000W/500W/ 600W/ 800W/1200W
Spot Size	: 12*13mm (300W/ 500W/ 600W) 12*23mm (800W/1200W) best size for light focusing
Crystal temperature	: ~20~3°C
Power requirement	: AC220V/110V 50~60Hz
NW	: 40 KG
GW	: 61 KG
Machine Size	: 40*35*110cm (L*W*H)
Package Size	: 57*52*121cm (L*W*H)



FEATURES:

1. 3 Wavelength : Treat all kind of hair color from black hair to white hair
2. Treat all skin types from white to dark skin
3. No pain and shorter treatment sessions
4. Effective and safe treatment for permanent hair removal



Naked Device



Injecting Gun (Handle)



Device Components Kit



- Deep clean
- Lifting tightening
- re-hydration
- freckles fading skin rejuvenation
- noninvasive skin beauty

SPECIFICATION:

Type	: Mesotherapy Gun
Certification	: E TUV ISO13485 FDA
Brand Name	: DFBEAUTY / DFLASER, DFBEAUTY
Model Number	: DF - MSG
Needle material	: Titanium Alloy
Treatment Head	: 3size micro needle
Injection depth	: 3.0mm
Technology	: RF mesotheroy



Before

After

Before

After



DERMA PEN WITH NEEDLE





Electromagnetic Puls

- Stimulates motor neurons to control pelvic floor muscles, to make pelvic floor muscles hyper-contract.
- Mild disorders: Vaginal wall relaxation, Genital dryness, Pain during sex; Difficulty in orgasm; Decreased sexual pleasure, Frequent urination, Constipation.
- In aggravated patients, Uterine prolapse, Prolapse of the front and rear walls of the vagina, Urinary incontinence.

The role of pelvic floor muscles in urinary incontinence

Pelvic floor muscles control the tissues of the urethra

Causes of Urinary Incontinence: Natural Aging, Childbirth, Menopause etc.



Effective treatment of urinary incontinence requires strengthening of pelvic floor muscle

Focused pulse of electromagnetic energy => deep penetration

=> stimulation of the entire pelvic floor area =>

one treatment=> 12,000 pelvic floor muscle contractions

Indications	: Stress urinary incontinence, urge urinary incontinence, mixed urinary incontinence
Course of treatment	: 6-8 times
Treatment interval	: 48 hours
Each treatment time	: 28 minutes
Duration of efficacy maintenance	: at least six months

KEGEL EXERCISER



A small comfortable probe is placed just a few centimetres into the vagina which sends gentle impulses to the correct pelvic floor muscles. This precise stimulation causes them to contract and get stronger quickly.

- **Strengthens pelvic floor after childbirth**
- **Helps and stress incontinence naturally**
- **Improve sexual intimacy & satisfaction**
- **Discreet pelvic floor toning therapy**
- **Flattens tummy**



INTIMATE AREA FILLER

with Lidocaine
WINK FEEL

**Second Diamond
to change your life**



Glycopeel

Glycolic Acid



Available strengths
20%, 35%, 50%, 70%

Indication :

- Aging skin • Stimulate collagen production • Skin rejuvenation
- Fine lines • Age spots
- Pigmentation • Acne

50ml

SALIC PEEL

Salicylic Acid



Available strength
30%

Indication :

- Salicylic acid exfoliates skin
- Improve skin texture and color.
- Oily Skin

50ml

MASTER PEEL

Master Peel exfoliating treatment



TCA 30% & Phytic colour (blue).

Indication :

- Pigmentation • Aging Skin
- Completely peels & removes dry flaky skin

100ml

Ablaze Peel



Glycolic acid 50 %
Lactic acid 10%
Salicylic acid 4%

- Restores a radiant complexion
- Acts as epidermis stimulator
- Regulates oily skin
- Dilated pores
- Improves fine lines and wrinkles
- Even skin tone

50ml

FUSION PEEL A
Lactic 45% + Citric 5% + Pyruvic 5% + Arginine 5%



50 ml

- Wrinkles • Dark circle • Aging skin

FUSION PEEL B
Lactic 15% + Kojic 5% + Citric 5% + Sal 2%



50 ml

- Skin lightening • Melasma
- Photodamaged skin

FUSION PEEL C
Sal 17% + Lactic 17% + Citric 8% + Gly 15%



50 ml

- Anti Acne • Anti Aging
- Freckles

FUSION PEEL D
Sal 20% + Man 10%



50 ml

- Stretch Marks • Scarring
- Post - inflammatory acne pigmentation

FUSION PEEL E
Gly 35% + Kojic 10% + Pyruvic 2% + Lactic 10%



50 ml

- Skin Rejuvenation • Superfacial
- Wrinkles

LA PEEL

ARGILAZ PEEL Lactic 20%, Arginine 20%



Lactic 20%
Arginine 20%

Applications:
Reduces Wrinkles
Increases Skin
Hydration level
Stimulate Collagen

50 ml

LA 40

LACTIC 40%



Lactic 40%

Applications:
Treats Hyperpigmentation,
Dehydrated skin,
Premature aging,
Skin Lightening,
Lightens Sun Damaged
skin, Fades age spots,
Acne.

50 ml

LA 80

LACTIC 80%



Lactic 80%

Applications:
Treating Dry, Rough &
Scaly skin, Excellent
Exfoliant and
Increases Skin
Hydration Level

50 ml

PEEL M1

MANDLAZ - 20 Mandelic 20%



Mandelic 20%

Applications:
Exfoliates Dead skin
cells, Reduces the
Apperance of
Pores Prevent
Breakout of Acne

50 ml

PEEL M2

MANDLAZ - 30 Mandelic 30%



Mandelic 30%

Applications:
Simulates Collagen,
Indicates for Oily and
Seborrheic skin,
Stimulates Collagen and
Proteoglycan Synthesis.

50 ml

PEEL M3


MANDLAZ - 40 Mandelic 40%



Mandelic 40%

Applications:
Fine lines, Crystic acne,
Large Pores, Darker
Skin Tones and Sun
Damaged skin,
Loosen Blackheads.

50 ml



Always in pursuit of innovative technologies pvt.ltd

Analyticals Medical Technologies Pvt.Ltd.

Shree Ambika Heritage, 1st Floor, Plot No.1, Sector 1,
Near Kharghar Railway Station, Navi Mumbai 410210

www.amtmedicals.com

info@amtmedicals.com

 +91 8976074090